

A photograph of the Nevada State Health Laboratory building. The building is a single-story, light-colored concrete structure with a flat roof. The words "NEVADA STATE HEALTH LABORATORY" are mounted on the upper portion of the facade in large, dark, sans-serif capital letters. Below the text, there is a series of windows and columns. The building is partially obscured by the green leaves of a tree in the foreground. The ground in front of the building is covered in green grass, with dappled sunlight and shadows from the tree. The sky is visible through the branches of the tree.

NEVADA STATE HEALTH LABORATORY

Foodborne Outbreak Investigation

Laboratory Component

Nevada State Health Laboratory

Overview

- Criteria for Testing
- Key Items for a Field Sample Kit
- Collection of Clinical Specimens
- Collection of Food Samples and Environmental Sampling
- Lab Requisition Forms and Chain of Custody
- Testing performed at the *NSHL* and Reference Laboratories
- *Nevada State Health Laboratory* Contact Information

Criteria for Testing

- Communication and team work is the key for a successful outcome of an investigation of a foodborne outbreak
- Environmental Health Assessment
- Epidemiological Investigation
- Laboratory Guidance and Support
- Education and Training
- Testing performed is determined by the local health department in collaboration with the above specialty areas as a foodborne outbreak investigation team

Criteria for Testing (Continued)

- Joint evaluation of the outbreak by the Response Team will direct what samples are to be collected and submitted to the lab for testing
- Information in the initial stages of an outbreak may be insufficient for an adequate assessment
- Accordingly the response team will make decisions of what to collect and then hold samples in the Laboratory until more information and assessment is available to direct testing
- Such a plan will ensure that the opportunity to collect the appropriate specimens is not missed while still ensuring that costly and time consuming lab resources are utilized effectively

Key Items for a Field Sample Kit

- Collection of Clinical Samples
- Collection of Food Samples
- Collection of Environmental Samples

- Be prepared

Collection Kit for Clinical Samples

- Cary-Blair Preservative Stool Sample Containers
- Sterile Specimen Containers without Preservative
- Ova & Parasite Collection Containers – Formalin and PVA
- Labels and marker pens
- *NSHL* Lab Requisition Forms
- Written Instructions for the proper collection and transport of clinical, usually stool, samples

Collection Kit for Food Samples

- Sterile Plastic Collection Containers
Such as Whirl Pac – Various sizes
- Sterile Disposable Collection Utensils
Spoons and Scoops
Tongs
Knife and cutting utensils
- Sterile Disposable Gloves
- Labels and Marker pens
- Chain of Custody Forms
- *NSHL* Requisition Forms

Collection Kit for Environmental Samples

- Sterile Collection Containers
Sterile Plastic such as Whirl Pac Bags
- Sterile Saline, Buffer or Enrichment Broth
- Sterile Swabs or Sponges
- Sterile Disposable Gloves
- Disposable Biohazardous or Garbage Bags
- Labels and Marker Pens
- *NSHL* Requisition Forms

Instructions for Collection and Transport of Clinical Samples

- Stool samples
- Emesis
- Completing lab requisition forms

Instructions for Collection of Stool Samples

- Do not pass directly into specimen containers
- Collect with newspaper on toilet set
- Add portion of stool sample to the appropriate specimens containers – do not fill above line
- Secure lid – No leaking samples
- Provide verbal and written instructions to patients
- Label specimen
Patient name, date and time of collection
Indicate if case patient or food handler
- Complete *NSHL* lab form
- Transport Cary-Blair specimen at room temperature without delay
- Transport sterile container without preservative at refrigerator temperature without delay (Norovirus)
- Transport sterile container without preservative frozen without delay (electron microscopy)
- Transport Ova & Parasite specimens at ambient temperature without delay
- Without delay is defined as delivery to the lab with in 24 hours of collection

Collection of Emesis

- Sterile container without preservative
- Label and transport refrigerated without delay
- Norovirus testing
- Not recommended for routine bacterial culture
- Requires immediate neutralization of specimen acidic pH at the time of collection

Instructions for Collection of Food and Environmental Samples

- Our results are only as good as the Samples we receive
- Contact the laboratory prior to collection and transport of samples for testing
- Ensures correct collection and transport and faster lab turn around times

Collection of Food Samples

- Wear sterile disposable gloves
- Use sterile technique
- Use sterile utensils for collection
- Use sterile transport containers
- Secure sample container – No leaking specimens
- Label all samples
 - Description and condition of food sample
 - Initials of sample collector
 - Time and date of collection

Collection of Food samples

Completing Lab Requisition Form

- Name of individual collecting sample
- Name of Local Health Department conducting investigation
- Contact information
- Date and time of collection
- Temperature, pH and any special conditions of the food at the time of collection
- Lot numbers, producer identification, site or origin of sample collection
- Signs, symptoms, onset times and number of persons ill
- Indicate Lab Tests to be performed

Collection of Food Samples

- Recommend minimum 250 grams of solid food
Prefer 500 grams
- Recommend 250 minimum milliliters of liquid food
Prefer 500 milliliters
- In general, samples should be collected from the geometrical center of the food item
- If in doubt, send the entire food sample
- Commercially packaged food samples must be retained and transported the original package
- Most food items are transported at refrigerator temperature
- If the food is already frozen at the time of collection, then transport the sample frozen

Collection and Transport of Food Samples

Sample	Transport
Solid Food	Sterile labeled container Refrigerate
Liquid Food or Beverage	Sterile labeled container Refrigerate
Frozen Food	Keep frozen - Label Send entire packaged sample
Canned Food	Label can or jar Transport at ambient temperature
Dry or Dehydrated Food	Label container Transport at ambient temperature
Meat or poultry	Place meat in sterile container Refrigerate

Instructions for Collection of Environmental samples

- Wear sterile gloves
- Use sterile technique
- Pre-wet sterile swab or sponge with sterile saline, buffer or enrichment broth
- Swab or sponge a one foot square environmental surface area
- Place swab or sponge in a sterile container with sterile saline, buffer or enrichment broth
- Label and transport without delay
- Complete *NSHL* Requisition Form

Lab Requisition Forms

- Complete requested information
- Sample collector and local health agency conducting investigation
- Date and time of collection
- Physical description and condition of sample
- Site of collection
- Lot numbers
- Signs, symptoms, on-set times, number of persons ill
- Indicate lab tests to be analyzed

Food-Borne Illness Requisition Form



NEVADA STATE HEALTH LABORATORY
 1660 NORTH VIRGINIA STREET
 RENO, NV 89503
 Phone: (775) 688-1335 Fax: (775) 688-1460
 MEDICAL DIRECTOR: L.D. BROWN, MD MPH

Client: _____ Clinic: Food Investigation
 Address: _____ Attn: _____

Phone: _____ Copy To: _____

Date Submitted: _____	Submitted By: _____
Date/Time Collected: _____	Phone Number: _____
Source: _____	Report To: _____
Complainant Name: _____	
Complainant Address: _____	
Complainant Symptoms (Check Box):	
<input type="checkbox"/> Fever (highest Temp. _____)	<input type="checkbox"/> Blurred Vision
<input type="checkbox"/> Cramps/Abdominal Discomfort	<input type="checkbox"/> Diarrhea
<input type="checkbox"/> Vomiting	<input type="checkbox"/> Nausea
<input type="checkbox"/> Headache	<input type="checkbox"/> Anorexia
<input type="checkbox"/> Severe Weakness	<input type="checkbox"/> Malaise
<input type="checkbox"/> Other (specify): _____	<input type="checkbox"/> Scleral Icterus (yellow eyes)
Date/Time Food Ingested: _____	Date/Time of Illness Onset: _____
Duration of Illness: _____	Number of People Ill: _____
Number of People Consuming Suspect Food (ill or not ill): _____	
Physical Description of Sample Being Submitted: _____	

**** Check box to the left of test to be ordered ****

BACTERIAL

<input type="checkbox"/> Bacillus cereus	<input type="checkbox"/> E. coli O157:H7	<input type="checkbox"/> Streptococcus group A
<input type="checkbox"/> Brucella species	<input type="checkbox"/> Listeria monocytogenes	<input type="checkbox"/> Vibrio species
<input type="checkbox"/> Campylobacter species	<input type="checkbox"/> Salmonella species	<input type="checkbox"/> Yersinia
<input type="checkbox"/> Clostridium botulinum	<input type="checkbox"/> Shigella species	
<input type="checkbox"/> Clostridium perfringens	<input type="checkbox"/> Staphylococcus aureus	

VIRAL

<input type="checkbox"/> Hepatitis A	<input type="checkbox"/> Norwalk Family of viruses	<input type="checkbox"/> Astrovirus, calcivirus, others
<input type="checkbox"/> Other (specify): _____		

PARASITE

<input type="checkbox"/> Cryptosporidium parvum	<input type="checkbox"/> Cyclospora cayetanensis	<input type="checkbox"/> Giardia lamblia
<input type="checkbox"/> Trichinella species	<input type="checkbox"/> Other (specify): _____	

MISCELLANEOUS

<input type="checkbox"/> Fungal (specify) _____	<input type="checkbox"/> Foreign Material Identification
---	--

Chain of Custody

- Legal document verifying integrity of the sample from the moment of collection in the field until the final deposition of the sample
- Official Sample
Integrity of the sample is verifiable
includes internal and external transfer information
- Unofficial Sample
Integrity to the sample is not verified
example – left-overs from an individuals refrigerator
- 3 most important things when conducting testing for a food borne outbreak
Document, Document and Document

Testing Performed at the *Nevada State Health Laboratory*

- Traditional Microbiology Cultures, Virology and Parasitology
- Toxin Assays
- Food Chemistry
- Molecular Biology Methods
- Reference Laboratory Testing

Foodborne Pathogens Tested for on Clinical and Food Samples at the *NSHL*

- *Salmonella* serotypes
- *Shigella* species
- *E.coli* O157:H7
- *Listeria monocytogenes*
- *Vibrio* species
- *Yersinia enterocolytica*
- *Staph aureus*
- *Streptococcus* group A
- *Bacillus cereus* (food samples only with enumeration)
- *Clostridium perfringens* (food samples only with enumeration)

Foodborne Pathogens Tested on Clinical Samples at the *NSHL*

- Norovirus – PCR
- Hepatitis A Virus - Serological
- Parasite – direct observation of parasite with special stains

Cyclospora cayatensis

Cryptosporidium parvum

Giardia lamblia

Entamoeba histolytica

Trichinella spiralis

Toxin Assays

- Staphylococcal Enterotoxin B
Test currently available – TRF method
- *Clostridium botulinum* Neurotoxin
Test validation in progress

Food Chemistry

- Aluminum, Antimony, Arsenic, Barium, Beryllium, Bismuth, Boron, Cadmium, Calcium, Cesium, Chloride, Chromium, Cobalt, Copper, Cyanide, Fluoride, Iron, Lead, Magnesium, Magnesium, Manganese, Mercury, Molybdenum, Nickel, Nitrate+Nitrite as N, Nitrite-N, Perchlorate, pH, Temperature, Potassium, Selenium, Silica, Silver, Sodium, Strontium, Sulfate, Thallium, Tin, Titanium, Total Phosphorus, Tungsten, Turbidity, Uranium ICP-MS, Vanadium, Zinc

Molecular Testing Currently Available at the *NSHL* PCR

- PCR based (Polymerase Chain Reaction) amplification of target Nucleic Acid

Norovirus – Clinical Specimens only

Salmonella – Food and Environmental samples

E.coli O157:H7 – Food and Environmental samples

Listeria species *and Listeria Monocytogens*

- Food and Environmental samples only

Molecular Testing Currently Available at the *NSHL*

PFGE

- Pulsed Field Gel Electrophoresis (PFGE)
Molecular sub typing of bacterial isolates for establishing laboratory links for epidemiology
- Currently available
Non-Typhoidal Salmonella
E.coli
Shigella
- Test validation in progress
Listeria monocytogens
Campylobacter
Vibrio species
- Each isolate PFGE pattern is posted on PulseNet and analysis performed for possible matches
- The *NSHL* participates in sending isolates to the CDC for NARMS

Reference Laboratory Testing

- The *Nevada State Health Laboratory* provides support for testing done at Reference Laboratories such as the CDC and FDA Labs
- Specimens and requisitions are submitted the *NSHL*
- Laboratory results are compiled and forwarded to the client and appropriate health agencies

Nevada State Health Laboratory

Contact Information

- 24 hour pager number 775-823-1150
- Fed Ex Account # 1046-0029-8
- Address
Nevada State Health Lab
1660 North Virginia Street
Reno, NV 89503
- Phone number 775-688-1335
- Fax number 775-688-1460



NEVADA
STATE HEALTH
LABORATORY

1000 W. WASHINGTON BLVD.